



Bromazepam – Ring Open or Ring Closed?

Jonas Malzacher, Lana Brockbals, Andrea Eva Steuer, Thomas Kraemer

Department of Forensic Pharmacology and Toxicology, Zurich Institute of Forensic Medicine, University of Zurich

When comparing the fragmentation of different analytes, two peaks were observed for bromazepam. While this can easily be explained by the typical instability of benzodiazepines, leading to ring-opened and ring-closed species, the simultaneous occurrence of these two species at the first of the two peaks deserved our attention. In-source recyclization was hypothesized as a likely explanation. Correspondingly, an increase in the ion source temperature resulted in an increase of the signal of ring-closed bromazepam at the retention time of the ring-opened species.




Fig.6: Extracted ion chromatograms of bromazepam and its ring-opened form, 20 hours after start of experiment (equilibrium was reached after approx. 13 h)




Fig.3: Reversible ring-cleavage reaction of bromazepam




Fig.4: Development of the area of ring-closed divided by ring-opened form




Fig.5: Development of the area of ring-closed divided by ring-open form

Methods

- Experiment conducted using HPLC-HRMS/MS approach
- Measurement employing Sciex ZenoTOF 7600 with reversed-phase chromatography
- Untargeted method in DDA mode
- Eluent A (10 mM ammonium formate:acetonitrile (99:1) + 0.1% formic acid) and eluent B (acetonitrile:ultrapure water (99:1) + 0.1% formic acid)
- Use of eluent A and eluent B with and without acid as solvents
- Analyte mix containing ten different benzodiazepines (bromazepam, pyrazolam, clonazepam, clonazolam, nitrazepam, diazepam, nordazepam, alprazolam, midazolam, triazolam) and two internal standards (lorazepam-d4 and trimipramine-d3)
- Sample measurement at ion source temperatures of 450 and 650 °C
- Repeated measurements in triplicates after incubation
- Higher density of measuring points at the outset
- Maximum incubation of 50 hours




Fig.1: Exemplary analytes with structural differences




Fig.2: Schematic representation of experimental set-up
Ten analytes, 100 ng/mL each (two deuterated analytes, 50ng/mL each)

Results

- Two bromazepam species (ring-opened and ring-closed) at the same RT due to alleged in-source recyclization
- Pyrazolam (also with bromo substituent) showed no ring opening
- Other evaluated benzodiazepines showed expected ring opening, however without subsequent in-source recyclization
- Signal intensity for bromazepam reduced by 50% at higher temperature (650 °C)
- Increased source temperature led to tenfold increase of ratio species 1 / species 2 (RT4.5)
- Equilibrium for bromazepam and corresponding opened form reached after approx. 13 hours
- Few analytes showed ring opening at higher pH, bromazepam not included

Conclusion

- Second peak for species 1 at 4.5 min associated to in-source cyclization reaction of the ring-open form, enhanced at higher ion source temperatures
- Bromine-substituted benzodiazepines with pyridine ring showed unique behavior
- Need for further investigation of other compounds of this class, especially those with similar structural elements
- Potential serious implications for quantification of benzodiazepines in biological samples (especially in the case of targeted quantification in scheduled MRM mode)

